Places to Start

Solvents

- Use brown borosilicate bottles to avoid algae growth
- Prepare solvent volume to be used up within 1 to 2 days
- Use only HPLC-grade solvents filtered through 0.2 µm filters

Preparing and powering up the pump

- Inspect solvent bottles and inlet filters for damage or coloring Always use seal wash when installed and purge the pump
- Use the appropriate system conditioning method

Daily tasks

- Replace aqueous and organic mobile phases every second day
- Check seal wash solvent
- Flush the system with the composition of your application

Weekly tasks

- Change seal wash solvent and bottle and inspect solvent filters
- Check system backpressure and change filters if necessary

Pump shutdown


- Flush all channels to remove salt deposits and particulate matter
- Flush the system with appropriate storage solvent and power down the system

Handling of acetonitrile

- If possible, use 5 to 10% of water in your mobile phase
- Be sure to avoid ACN evaporation
- Don't leave ACN on the system for more than 2 to 3 days – Perform a periodic warm water wash (60 to 70 °C) if you face problems







Maintenance

Agilent Lab Advisor software helps you manage your Agilent LC instruments to achieve high-quality chromatographic results in the most efficient way by ensuring high instrument performance, productivity, and reliability. It is available free-of-charge.

- Diagnostic tests to evaluate performance - Easier maintenance of all Agilent LC modules - Comprehensive reports generated to ease communication

- with Agilent service

Training courses are available at: https://www.agilent.com/crosslab/university

Discover more best practices for using an Agilent LC system https://www.agilent.com/chem/lc-best-practices

Retention Time Drift

Pressure Fluctuation

MM

Pressure Increase

P/

High Column

Backpressure

 P/ Λ

Drifting Baseline

Possible Cause	Solution
Inconsistent online mobile	Ensure gradient system delivers
phase mixing	constant composition; compare with manual preparation of mobile phase
Variation in column temperature	Thermostat or insulate column; ensure constant lab temperature
Insufficient equilibration time	Make sure at least 10 column
with gradient run or change in isocratic mobile phase	volumes pass through column after sample run
Selective evaporation of	Less vigorous helium sparging; keep
mobile phase component	solvent reservoirs covered; prepare fresh mobile phase
Contamination buildup	Occasionally flush column with strong solvent
Column overloaded	Decrease injection volume
with sample	or concentration
Possible Cause	Solution
Leak in the system	Identify the channel and clean
	or replace check valve; replace
	pump seals
Buildup of particulates	Filter sample and mobile phase
Bubble in pump	Perform solvent degassing; sparge solvent with helium
Possible Cause	Solution
System blockage	Check flowpath (needle seat,
	capillaries, filter and frits)
Water/organic systems:	Test buffer-organic mixtures
buffer precipitation	to ensure compatibility
Possible Cause	Solution
Column blockage	Better sample cleanup; use guard column
Mobile phase viscosity	Use lower viscosity solvents or
too high	higher temperature
Particle size too small	Use larger d_p packing
Plugged inlet frit	Replace column
Possible Cause	Solution
Positive/negative direction:	Flush column; clean up sample;
contaminant buildup/elution	use pure solvents
Positive/negative: difference in refractive index of injection solvent	Use mobile phase for sample solvent

Temperature changes



Insulate and thermostat column

and tubing






Get answers. Share insights. Join the Agilent Community at https://community.agilent.com

Use degassed HPLC-grade solvents;

Solution

with nase



Noisy Baseline


. / ^	Contamination	flush system; clean up sample
mm hour hour	Detector problems	Check number of hours of UV lamp; replace UV lamp or flow cell
Ghost Peaks	Possible Cause	Solution
A A A.	Peaks from previous injection	Flush column to remove contaminants; check with blank injection
	Contamination; unknown interferences in samples	Proper sample cleanup
	Ion pair: disequilibrium	Prepare sample in actual mobile phase to minimize disturbance
	Contaminated mobile phase	Check your mobile phase
	Bubbles in solvent	Check and degas your solvents
Peak Tailing	Possible Cause	Solution
$\bigwedge \bigwedge \bigwedge \to \bigwedge \bigwedge \bigwedge$	Unswept dead volumes	Minimize number of connections; ensure injector seal is tight; ensure fittings are properly seated
JUULJUUL	Column performance	Change mobile phase; replace column
	Silica-based: column degradation	Use specialty, polymeric, or sterically protected column
	Silica-based: basic interactions with stationary phase	Use stronger mobile phase or add appropriate base (e.g., TEA)
Peak Broadening	Possible Cause	Solution
A A A	Injection volume too large	Decrease injection volume or solvent strength of injection solvent; use gradient methods

Possible Cause

Contamination



Sensitivity Problems		



Leaks

Low sampling rate of data system Detector cell volume too large

Injection volume too large

Possible Cause Peaks are outside of sensitivity range of detector

Sample-related losses during preparation

Possible Cause White powder at fitting/ loose fitting System leak

olume or solvent solvent; use gradient methods Increase data rate Use smallest possible cell volume

Decrease injection volume

Solution Dilute/concentrate sample to bring into linear region

Use internal standard during sample preparation; optimize sample preparation method

Solution Tighten fittings; replace capillaries

Identify location checking leak sensors/errors; check flow cell

For Lab Advisor software, please visit: https://www.agilent.com/chem/lab-advisor





© Agilent Technologies, Inc. 2019 Printed in the USA, March 1, 2019 5994-0709EN