



### Kurzgewindeverschlüsse – 9mm PP bunt inkl. Septum

	PTFE virginal	Naturkautschuk/ TEF transparent	RedRubber/ PTFE beige	Silicon weiß/ PTFE rot	PTFE rot/ Silicon weiß/ PTFE rot	Silicon weiß/ PTFE blau	Silicon weiß/ PTFE rot
	53° shore D 0,2mm	60° shore A 1,0mm	45° shore A 1,0mm	55° shore A 1,0mm	45° shore A, 1,0mm	55° shore A 1,0mm geschlitzt	55° shore A 1,0mm angeschlitzt(Y)
 <b>Rosa Glitzer</b>	69040-0915	69041-0915	69042-0915	69043-0915	69044-0915	69045-0915	#
 <b>Transparent</b>	60981-0915	60478-0915	62011-0915	60481-0915	60480-0915	60852-0915	62021-0915
 <b>Blau</b>	60982-0915	60867-0915	61819-0915	60838-0915	60868-0915	60869-0915	62016-0915
 <b>Rot</b>	61337-0915	61176-0915	62012-0915	61178-0915	61177-0915	61179-0915	#
 <b>Schwarz</b>	61668-0915	61570-0915	62013-0915	61572-0915	61571-0915	61669-0915	#
 <b>Grün</b>	61539-0915	61356-0915	61911-0915	61332-0915	61485-0915	61746-0915	#
 <b>Gelb</b>	#	61542-0915	62014-0915	61527-0915	61486-0915	61745-0915	#
 <b>Pink</b>	#	#	62470-0915	62473-0915	#	62472-0915	#

VE = 1.000 Stück

Weitere Flaschen & Verschlüsse finden Sie in unserem Vials & Caps Katalog [>öffnen<](#)

DURATEC Analysetechnik GmbH | Rheinauer Straße 4 | D-68766 Hockenheim | [support@duratec.de](mailto:support@duratec.de) | [www.duratec.de](http://www.duratec.de)