



Kurzgewindeverschlüsse – 9mm PP bunt inkl. Septum

	PTFE virginal 53° shore D 0,2mm 	Naturkautschuk/ TEF transparent 60° shore A 1,0mm 	RedRubber/ PTFE beige 45° shore A 1,0mm 	Silicon weiß/ PTFE rot 55° shore A 1,0mm 	PTFE rot/ Silicon weiß/ PTFE rot 45° shore A, 1,0mm 	Silicon weiß/ PTFE blau 55° shore A 1,0mm geschlitzt 	Silicon weiß/ PTFE rot 55° shore A 1,0mm angeschlitzt(Y) 
 Rosa Glitzer	69040-0915	69041-0915	69042-0915	69043-0915	69044-0915	69045-0915	#
 Transparent	60981-0915	60478-0915	62011-0915	60481-0915	60480-0915	60852-0915	62021-0915
 Blau	60982-0915	60867-0915	61819-0915	60838-0915	60868-0915	60869-0915	62016-0915
 Rot	61337-0915	61176-0915	62012-0915	61178-0915	61177-0915	61179-0915	#
 Schwarz	61668-0915	61570-0915	62013-0915	61572-0915	61571-0915	61669-0915	#
 Grün	61539-0915	61356-0915	61911-0915	61332-0915	61485-0915	61746-0915	#
 Gelb	#	61542-0915	62014-0915	61527-0915	61486-0915	61745-0915	#

VE = 1.000 Stück

Weitere Flaschen & Verschlüsse finden Sie in unserem Vials & Caps Katalog [>öffnen<](#)

DURATEC Analystechnik GmbH | Rheinauer Straße 4 | D-68766 Hockenheim | support@duratec.de | www.duratec.de